

# TD-1241

## Non-Contact Forehead Thermometer



### FEATURES

- 3 Measuring mode: Forehead/Surface/Room temperature
- Accurate and Fast
- LCD Backlight



### SPECIFICATIONS

Dimension	150 (L) x 48.5 (W) x 55 (H)mm
Weight	125.8 g (With Battery)
Resolution	0.1°C/0.1°F
Measurement Range	Forehead: 22°C-44°C
	Surface: 0°C -100°C
	Room: 10°C -40°C
Accuracy	Forehead: 36°C-39°C±0.2°C(22°C-35.9°C, 39.1°C -44°C±0.3°C)
	Surface/Room: ±.1°C
Communication	No Data Transmission (TD-1241A) / Bluetooth smart (TD-1241B)
Measurement Units	°C(preset)/ °F
Power source	2*1.5V AA
Battery life	> 3000 measurements
Memory	30 measurements
Power saving	30 seconds
Operating condition	10°C to 40°C (50°F to 104°F), 95% RH or less (non-condensing)
Storage/Transportation condition	-20°C to 60°C (-4°F to 140°F), 95% RH or less (non-condensing)
Measuring Distance	3-7cm
Units per Carton	30
Net Weight per Carton	7.5kg
Gross Weight per Carton	8.5kg
Carton Size	42x39x41cm



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No.127, Wugong 2nd Rd., Wugu Dist.,  
New Taipei City 24888, Taiwan  
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288