

TD-8255

Fingertip Pulse Oximeter

Mb^H
Mobile Health



FEATURES



Large LCD
with backlight



Simple
Operation



Automatically
power-off function



URIGHT SpO2 Manager App



SPECIFICATIONS

| | |
|-------------------------|------------------------------------|
| Dimension | 63 (L) x 35 (W) x 32 (H) mm |
| Weight | 40g (without battery) |
| Heart Rate Range | 30 ~ 250 beats per minutes (bpm) |
| SpO2 Range | 0% ~ 100% |
| Accuracy (SpO2) | 80%~100% ±2 digits |
| Method | Dual wavelength LED |
| Accuracy (Heart rate) | ±1 bpm or ±1% whichever is greater |
| Power Source | 2 x AAA |
| Time to power saving | 15 seconds after removing finger |
| Backlight | LCD backlight |
| Communication | Bluetooth smart |
| Operating Temperature | 10°C ~ 40°C |
| EMC Standard | IEC60601-1-2 |
| Safety Standard | IEC60601-1 |
| Performance Standard | ISO80601-2-61 |
| Unit per Carton | 80 |
| Net Weight per Carton | 11.2kgs |
| Gross Weight per Carton | 12.5kgs |
| Carton Size | 518 (L) x 423 (W) x 300 (H) mm |



www.taidoc.com

TaiDoc Technology Corp.

B1-7F., No. 127, Wugong 2nd Rd., Wugu Dist.,
New Taipei City 24888, Taiwan
Tel: +886-2-6625-8188 Fax: +886-2-6625-0288